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**Internet: [www.tka.de](http://www.tka.de)**

**E - Mail: [info@tka.de](mailto:info@tka.de)**

### Table of dimensions and performance:

	TKA DI 150	TKA DI 150-DEN	TKA DI 300
<b>System including conductivity meter</b>			
Flow rate max.	50 l/h	50 l/h	50 l/h
Capacity at 10° TDS*	150 litres	2x 150 Liter	300 litres
Operating pressure	8,5 bar	8,5 bar	8,5 bar
Pure water quality	0,1 – 20 µS/cm	0,1 – 20 µS/cm	0,1 – 20 µS/cm
Water temperature max.	30°C	30°C	30°C
Electrical supply	230/50-60 V/Hz	230/50-60 V/Hz	230/50-60 V/Hz
Diameter x Height**	170 x 430 mm	2x 170 x 430 mm	170 x 705 mm
Housing	PP	PP	PP
Weight	3 kg	4,3 kg	4 kg
<b>Artikel No.</b>	<b>01.0150</b>	<b>01.0150-DEN</b>	<b>01.0300</b>
Consumable material required:			
<b>Cartridge filled with disposable resin</b>	2 pieces	2 pieces	2 pieces
<b>Artikel No.</b>	<b>01.0151</b>	<b>01.0151</b>	<b>01.0351</b>

\* Total dissolved solids

\*\* including the conductivitymeter

### Equipment supplied:

- \* Pressure resistant plastic housing with 3/4" connections
- \* Conductivity meter, analog, display 0 – 50µS
- \* Feedwater hose incl. flow restrictor; 1,5 m  
**(If you need to replace the feedwater hose you are obligated to order the original TKA hose with the flow restrictor.)**
- \* Pure water hose; 1,5 m
- \* Wall mounting
- \* Leather housing wrench for the opening of filter housing.

### Safety Precautions:

01. Using the cartridge in any other way than that described in these operating instructions invalidates the guarantee!.
02. Deionized water must on no account be consumed in any other way.  
**It is not fit for drinking.**

### Installation:

Protect the room in which the system is installed from frost. Minimum temperature + 2°C.

### On-site requirements:

1. An R 3/4" feedwater (rawwater) connection.
2. A DN 50 (38.5 mm internal diameter) floor drainpipe must be available.
3. There must be no source of heat in the vicinity of the system.
4. Connect only the hoses supplied of pure water resistant materials to the down-stream side of the plastic housing.

When nonferrous metal material is used, there is a possibility of decomposition over time.

### Insert of the cartridge:

1. Screw the bottom of the housing off.
2. Insert the cartridge with the gasket upwards.
3. Screw the bottom of the housing back on.

### Connecting the system:

1. Close the feedwater tap.
2. Connect the inlet hose with the straight connector (**flow restrictor integrated**) to the **water tap (the connector is labelled "rawwater")** and the angled connector to the **plastic housing connector marked „IN“**.  
Take care that the house is not buckled.
3. Screw the pure water hose onto the plastic housing connector marked **OUT**.
4. Connect the open end of the pure water hose to your user.
5. When used to supply sensitive applications, it may be necessary to install a filter to protect the user against particles from the ion exchanger.

**Important: Be sure not to mix up the feedwater and pure water hoses.**

### Vent the cartridge:

1. Loosen the house on the pure water outlet.
2. Open the feedwater tap.
3. As soon as water emerges, re-tighten the hose on the pure water outlet.
4. The system is now vented and ready for use.

The disposable cartridge must be replaced as soon as the limiting conductivity which you have specified is reached (reading on the conductivity meter).

## Changing the cartridge:

The cartridge is exhausted when the pointer of the conductivity meter leaves the green field, or shows 20  $\mu\text{S}/\text{cm}$ .

1. Turn off the supply of water. Take out the mains plug. Disconnect the feedwater hose from the water tap; allow water to run out of the hose.
2. Screw off the bottom of the housing.
3. Remove the exhausted cartridge from the bottom of the housing and dispose of it with normal waste (dustbin).
4. Insert the new cartridge with the gasket upwards.
5. Screw the bottom of the housing back on.
6. Re-connect the feedwater hose to the feedwater tap.
7. Vent the cartridge.
8. The system is now ready to use.

Distributed by:

**Carl Stuart Limited**

ADVANCED APPLIED TECHNOLOGIES

**Contact Us:**


Irl Ph: 01 4523432

UK Ph: 08452 30 40 30

Web: [www.carlstuart.com](http://www.carlstuart.com)

Email: [info@carlstuart.com](mailto:info@carlstuart.com)


## TKA DI 150

- |  |  |
|--|--|
| 1) Leitfähigkeitsmessgerät 0-50 $\mu$ S (analog)<br>Analog conductivity meter 0-50 $\mu$ S | Art.-Nr.: 26.0002<br>Art.-No.: 26.0002 |
| 2) Filtergehäuse 10", komplett<br>Plastic housing, complete                                | Art.-Nr.: 06.5310<br>Art.-No.: 06.5310 |
| 3) Ionenaustauscher-Einsatz, Set (2 Stück)<br>Ion exchange cartridge, 2 pieces             | Art.-Nr.: 01.0151<br>Art.-No.: 01.0151 |
| 4) Rohwasserschlauch inkl. Dosierblende<br>Feedwater hose, incl. flow restrictor           | Art.-Nr.: 25.0036<br>Art.-No.: 25.0036 |
| 5) Reinwasserschlauch<br>Purewater hose  | Art.-Nr.: 18.0042<br>Art.-No.: 18.0042 |


## TKA DI 150-DEN

- |    |   |  |
|----|---|--|
| 1) | Leitfähigkeitsmessgerät 0-50 $\mu$ S (analog)<br>Analog conductivity meter 0-50 $\mu$ S | Art.-Nr.: 26.0002<br>Art.-No.: 26.0002 |
| 2) | Filtergehäuse 10", komplett<br>Plastic housing, complete                                | Art.-Nr.: 06.5310<br>Art.-No.: 06.5310 |
| 3) | Ionen-austauscher-Einsatz, Set (2 Stück)<br>Ion exchange cartridge, 2 pieces            | Art.-Nr.: 01.0151<br>Art.-No.: 01.0151 |
| 4) | Rohwasserschlauch inkl. Dosierblende<br>Feedwater hose, incl. flow restrictor           | Art.-Nr.: 25.0036<br>Art.-No.: 25.0036 |
| 5) | Reinwasserschlauch<br>Purewater hose  | Art.-Nr.: 18.0042<br>Art.-No.: 18.0042 |

TKA DI 300

- |    |   |  |
|----|---|--|
| 1) | Leitfähigkeitsmessgerät 0-50µS (analog)<br>Analog conductivity meter 0-50µS   | Art.-Nr.: 26.0002<br>Art.-No.: 26.0002 |
| 2) | Filtergehäuse 20", komplett<br>Plastic housing, complete                      | Art.-Nr.: 06.5320<br>Art.-No.: 06.5320 |
| 3) | Ionenaustauscher-Einsatz, Set (2 Stück)<br>Ion exchange cartridge, 2 pieces   | Art.-Nr.: 01.0351<br>Art.-No.: 01.0351 |
| 4) | Rohwasserschlauch inkl. Dosierblende<br>Feedwater hose, incl. flow restrictor | Art.-Nr.: 25.0036<br>Art.-No.: 25.0036 |
| 5) | Reinwasserschlauch<br>Purewater hose  | Art.-Nr.: 18.0042<br>Art.-No.: 18.0042 |