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**Read these operating instructions carefully before installing
and starting-up the system!**

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ECC-Declaration of Conformity

**acc. to EEC-Directive 98/37/EEC
- Machine Directive -**

We hereby declare that the products named below conform in their design and construction, as well as in the versions we have brought to market, to the fundamental safety and health requirements of **EEC-Directive 98/37/EEC**.


This declaration is invalidated should changes not agreed to by us are made to the machine.

Designation of the machine: **Ion exchange System**
Machine type: **TKA DI 1500 – TKA DI 15000**
Article number: **02.1500 – 02.15000**
Applicable EEC Directives: **EEC Machine Directive (98/37/EEC)**
Standards applied: **DIN EN ISO 12100-1, DIN EN ISO 12100-2, DIN EN 1050**

Niederelbert, 12. May 2003

Authorized representative of the manufacturer

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- 1 Feedwater connector, male R $\frac{3}{4}$ "
 - 2 Feedwater hose, R $\frac{3}{4}$ " female nut, length 1.5 m
 - 3 Deionized water hose, R $\frac{3}{4}$ " female nut, length 1.5 m
 - 4 Measuring cell and T-distribution
 - 5 TKA Control 330 conductivity meter incl. wall mounting
 - 6 Transformer plug
- Item 2 - item 6 = extend of delivery **Art.-No: 04.1805**
- 7 Solenoid valve, chem. nickel plated, in Art-No. **04.1807** included
 - 8 TKA Control 330 conductivity meter with switch socket for connection of a solenoid valve **Art.-Nr.: 04.1807**
 - 9 TKA Control 330 conductivity meter with potential-free-contact **Art.-Nr.: 04.1806**
- 10 TKA conductivity meter, analog
 - 11 Gasket $\varnothing 24/\varnothing 12 \times 2$
- Pos. 2, 3, 10, 11 = extent of delivery **Art.-No: 04.1601**

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General information

1. Read these operating instructions carefully before installing the cartridge system and putting it to operation.
2. The manufacturer will accept no liability when the system is operated improperly, or for other than the intended purpose.
3. The CE-mark ceases to be valid when constructional changes are made to the system, or products from other manufacturers incorporated.
4. Protect the system from frost. The temperature in the room in which the system is installed must be at least +2°C.
5. The installation area must have a drain at floor level with at least DN 50 pipe (40 mm o.d.). The manufacturer is not liable for damage caused by overflow of water. Should no such drain be available, we recommend that a water watcher (article no.: 16.0129) be installed.
6. The water supply must be turned off when the system is out of use for any length of time (over night, weekends, holidays) to avoid leakage.
7. The water demineralization cartridge must only be fed with cold water. The sealing material is not covered by the guarantee at temperatures above 32° C.
8. Follow the laws and regulations, including the accident prevention regulations, which are in force for the location where the system is to be installed.
9. When used to supply sensitive applications, it may be necessary to install a filter to protect the user against particles from the ion exchanger.
10. The guarantee is valid for 12 months.

For your security: take care of this general information!

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Important:

**The mixed-bed water demineralizers
TKA-DI 1500 - DI 15000 are
solely for the
purpose of purifying
drinking water (water which
fulfils the regulations for drinking water).**

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Purified water is required for many procedures and production processes.

TKA mixed-bed water demineralizers produce pure, demineralized water by means of ion exchange. They remove all dissolved salts, as well as carbonic and silicic acids, from the water.

Particularly frequent users:

- Hospitals
- Clinics and Sanatoriums
- Pharmacies
- Laboratories
- Industry

Application areas:

- Steam generation for sterilizers
- Steam generation for heating and for ironing machines
- Rinse water for laboratory washers and glassware rinsers
- In the electronic industry for the production of components
- Water for electroplating
- In laboratories and health care
- Air humidifiers and air conditioners
- For re-filling accumulators
- As polishing filters following reverse osmosis systems

Description of the function:

TKA mixed-bed water demineralizers contain a carefully chosen mixed-bed resin, made up from a strong cation and a strong anion ion exchange resin.

The fine mixing of these results in a very high quality water treatment unit.

During use, the feed water (drinking water) flows through the ion exchange resins from the top downwards. The water passes by innumerable chains of the resins, whereby salts are removed. The degree to which the water is desalted is shown by a conductivity meter in the purified water piping. When the conductivity limit which has been set down for the particular application is reached, the mixed-bed water demineralizer must be replaced and regenerated.

Parts supplied:

Mixed-bed water demineralizer consisting of an AISI 316 stainless steel housing filled with ion exchange resin. The housing has R3/4" stainless steel connectors and is

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Determination of the exchange capacity

With a feedwater quality of 1°dH = 1.25°eH total dissolved solids (approx. 30 µS), the exchange capacity is 10 times the number given in the catalogue number. Divide this value by your actual TDS content and you have the exchange capacity for the chosen TKA system.

Example: Type DI 2800: with water of 1°dH TDS = exchange capacity of 28,000 liters
 with water of 10°dH TDS = exchange capacity of 2,800 liters

Note: A high carbonic acid and silicat content in the water causes a decrease in the exchange capacity.

Accessories

- **A set of connecting hoses**
with female R ¾" spouts, PVC
Catalogue no: 04.1690
- **A set of connecting hoses**
with quick connect coupling and R ¾" adapter
Catalogue no.: 04.1691
- **Conductivity meter analogous**
to be installed into a front panel
Catalogue no.: 04.1600
- **Conductivity meter analogous**
to be mounted onto the cartridge
Catalogue no.: 04.1601
- **Conductivity meter digital**
Conductivity meter digital Control 330, with limit value switching up to 30 µS/cm and acoustic warning and optical signal.
Catalogue no.: 04.1805
Conductivity meter digital Control 330, with limit value switching up to 30 µS/cm, acoustic warning and optical signal, and potential-free contact.
Catalogue no.: 04.1806
Conductivity meter digital Control 330, with limit value switching up to 30 µS/cm, acoustic warning and optical signal, solenoid valve and potential-free contact.
Catalogue no.: 04.1807

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Installation

The room in which the system is installed must be protected from frost and must have a minimum temperature of +2°C.

Requirements at the site are as follows:

1. A feedwater connection R3/4"
2. A floor drain with DN 50 (38.5 mm internal diameter) piping.
3. There must be no source of heat in the vicinity of the area where the unit is to be installed.
4. Only the hoses supplied, or material which is compatible with purified water, are to be used downstream of the water demineralizing cartridges.

When nonferrous metal is used, subsequent corrosion cannot be ruled out.

Connection:

- Close the feedwater supply tap.
- Remove the caps from both connectors. (Keep the caps and the packaging, you will need it when returning the cartridge for regeneration).
- Without exception, lay the flow restrictor gasket for your type of water deionizer in the socket of the hose that is to be connected to the feedwater tap. (Applies only to types DI 1500 to DI 4000)

■ **Important! The insertion of the flow restrictor gasket is an absolute necessity for proper functioning of the water deionizing cartridge!**

- Tighten the spout with the flow restrictor inside onto the closed feedwater supply tap.
- Tighten the other end of the feedwater hose onto the connector of the water demineralization cartridge which is marked **EINGANG (INLET)**.
- Tighten the purified water hose onto the connector of the water demineralization cartridge which is marked **AUSGANG (OUTLET)**.
- Connect the other end of the purified water hose via a conductivity meter to the user.

■ **Important: Do not mix up the feedwater and purified water hoses.**

- Information on the connection to a conductivity meter is on the front page of this Operating Instructions.

Venting the cartridge

- Open the vent valve in the central closing lid of the cartridge by turning it anti-clockwise.
- Open the feedwater supply tap and wait until water comes out of the valve.
- As soon as water appears, close the valve and the feedwater supply tap.
- The system is now ready for use.

When the conductivity limit which you have decided upon has been reached (shown by the conductivity meter), the mixed-bed water demineralizer must be regenerated.

Regeneration

In order that you always have a supply of purified water, we suggest that you purchase an additional cartridge. When the cartridge which is exhausted is removed and sent for regeneration, it can be simply replaced by the other cartridge.

Please send exhausted cartridges to the following address:

TKA Water Purification Systems GmbH
Stockland 3
D-56412 Niederelbert
Germany

Telephone: +49 (0) 26 02 1 06 99-0

Fax: +49 (0) 26 02 1 06 99-50

www.tka.de

eMail: info@tka.de

Specifications given are not binding and we preserve the right to change them.

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Carl Stuart Limited
TKA Wasseraufbereitungssysteme GmbH