Distributed by:

**Lab Unlimited**  
CARL STUART GROUP

Tallaght Business Park  
Whitestown, Dublin 24,  
Ireland  
D24 RFK3

Tel: (01) 4523432

Fax: (01) 4523967

E-mail: [info@labunlimited.com](mailto:info@labunlimited.com)

Web: [www.labunlimited.com](http://www.labunlimited.com)


Quatro House, Frimley Road,  
Camberley,  
United Kingdom  
GU16 7ER

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Dear customer,  
Congratulations on the purchase of an  
ebro product.

We hope than you can make use of this  
product for a long time and that it helps  
you with your work.

Information that are useful and important  
for understanding the function are high-  
lighted in the instructions text.

**Standards:**

 The conformity certificate con-  
firms that this product fulfills  
CE guidelines.

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## Contents

<b>Overview .....</b>	<b>7</b>
<b>Safety information .....</b>	<b>9</b>
<b>What is the meaning of the term ...</b>	<b>11</b>
<b>Unpacking / scope of delivery .....</b>	<b>13</b>
Available accessories .....	13
<b>Turning the device on / of .....</b>	<b>15</b>
<b>Taking measurements .....</b>	<b>17</b>
Setting the temperature .....	17
Carrying out the measuring process ..	17
<b>Storing the currently measured value .</b>	<b>19</b>
<b>Changing device settings .....</b>	<b>21</b>
<b>Changing the sensor .....</b>	<b>25</b>
<b>Calibrating the device .....</b>	<b>27</b>
Preparing for calibration .....	27
Calibration .....	29
<b>What to do when ... ..</b>	<b>35</b>
Changing the battery .....	35
<b>Cleaning and care .....</b>	<b>38</b>
<b>Disposal .....</b>	<b>39</b>
<b>Technical data .....</b>	<b>41</b>
<b>Conformity declaration .....</b>	<b>42</b>

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
## Overview

The PHT 810 can be used to measure the pH value of a liquid.

Various sensors are available for different applications.

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## Safety information

**Do not use the device in explosive areas. Danger to life!**

**Do not bring the device into an environment warmer than 60°C. Lithium battery explosion risk!**

**In order to avoid damage to the device and the sensor, please also observe the following:**

- Protect the device **housing** from direct contact with water.
- Always install the sensor according to the instructions.

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## What is the meaning of the term...

### ... BNC connector

A BNC connector is a coaxial connector with a so-called bayonet nut. Coaxial means that an inner conductor is surrounded by a hose-shaped outer conductor, which has a shielding effect due to its shape. A bayonet nut operates according to the principle „put on - turn slightly (usually approx. 90°) – connected“.

### ... pH-value

The pH value is a measurement for the concentration of hydrogen ions in a solution. The value is stated as a negative base 10 logarithm ( $=10^{-\text{pH}}$ ). Acidic solutions have pH values  $< 7$ ; alkaline solutions have pH values  $> 7$ . „pH“ comes from „**p**ondus **H**ydrogenii (lat.: pondus = weight, hydrogenium = hydrogen).

### ... Buffer solution

A buffer solution absorbs contaminants with acids or bases without a material change in the pH value. Therefore the pH value of a buffer solution ideally remains constant, even when an acid or base is added (e.g. through immersion of a sensor coated with an acid or base).

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## Unpacking/Scope of delivery

Check the packaging and contents for sound condition upon receipt. Also verify that the package contents match your order.

In addition to the device with sensor and this user manual, various accessories can be included in the delivery:

### **Available accessories:**

Description	Identification
Case for PHT 810	AT100
Electrode for measurements in water, galvanic..	AT200
Electrode for measurements in chemistry, pharmaceuticals	AT201
Electrode for cut-in measurements in foods	AT202
Electrode for cut-in measurements in foods, with plastic shaft	AT204
Electrode for cut-in measurements in foods, with plastic shaft and grip	AT205
Punch	AT300
Buffer solution pH4 (50ml)	AT400
Buffer solution pH7 (50ml)	AT401
KCL solution (50ml)	AT405
Electrode cleaner	AT410
Drop protection	AG140

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**Should you have reason for a complaint, please contact us. Our contact information is found on the back of the envelope.**

## Turning the device on/off

◆ To turn the device on, press the »ON/HOLD« key.

After a few seconds (self-test) the device is ready to take measurements and already displays the current pH value.

A small battery symbol that shows the remaining charge in the device battery (see page 32) appears at the top left of the display.

◆ To turn the device off, press the »ON/HOLD« key until „OFF“ appears on the display (several seconds).

If you forget to turn off the device it will turn off automatically – provided this setting has not been modified (see page 18) – after two hours of operation to conserve the battery.

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
## Taking measurements

In order to achieve the highest level of measurement accuracy, you need to determine the temperature of the substance to be measured and set this on the device each time before you measure the pH value.

### Setting the temperature

- ◆ Briefly press the »CAL« key.
- ◆ Press the keys ▼ or ▲ until the actual temperature of the substance to be measured is shown on the display.
- ◆ Press the »ON/HOLD« key.

### Carrying out the measuring process

- ◆ Insert the sensor into the substance to be measured and read the displayed pH value.

**If the display shows the message »HI« or »LO«, please see page 32, „What to do when...“.**

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
## Storing the currently measured value

In order to better recognize changes in the measured value, you can store the current value on the display.

◆ In order to store the currently measured value on the display, briefly press the »ON/HOLD« key one.

A second display line now appears underneath the currently measured value. This contains the label HOLD and, to the right of this, the stored pH value.

The measured value in the upper display line is still continuously refreshed (see page 20, „Rate“, for the refresh rate).

The stored pH value is retained until you press the »ON/HOLD« key again.

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## Changing the device settings

You can set:

- Whether the device displays the temperature in °Celsius or °Fahrenheit.
- Whether the setting to automatically power down the device is activated or not.
- How often the measured value is refreshed (1 to 15 seconds).

In order to change device settings you have to call up the «User Menu». Proceed as follows:

### **Assumption: The device is turned off.**

- ◆ Press the »ON/HOLD« key until »PHT« is shown on the display.
- ◆ Now press the »CAL« key until »USER« is shown on the display.

The device is now in configuration mode and shows the «User Menu». This consists of the menu items »UNIT«, »SHDWN« and »RATE«.

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UNIT Temperature unit, °C or °F;

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SHDWN automatic shutdown after 2 hours of operation

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ON = shutdown activated  
OFF = shutdown deactivated

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RATE Refresh rate in seconds; can be set from 1-15 seconds; above 15 the scale reverts to 1.

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**Modified values are automatically stored right away.**

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◆ To leave the «User Menu» turn the device off by pressing the »ON/HOLD« key.

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## Changing the sensor

The sensor is connected to the device with a BNC connector. To change the sensor, proceed as follows:

- ◆ Turn the sensor on the BNC connector by 90° counter-clockwise while slightly pressing on the sensor in the direction of the device.

- ◆ Pull the sensor on the connector out downwards.

**Do not pull on the sensor!**

- ◆ Insert a new sensor.

- ◆ Turn the sensor and/or the extension cable on the BNC connector 90° clockwise.

**The BNC connector must lock in place!**

- ◆ Ensure the sensor connection is firmly seated.

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## Calibrating the device

In order to maintain a high level of measurement accuracy, the device must be calibrated every time the sensor is changed and after every three months of use at the latest. During calibration the gradient of the device is compared to one, two, or three buffer solutions that have a defined pH value.

**The larger the number of buffer solutions used and therefore the more points on the gradient are calibrated, the higher the subsequent measurement accuracy of the device will be.**

The device has a special calibration mode that guides you through the calibration.

### Preparing for calibration

To calibrate the device, you require at least one and at the most three buffer solutions with different pH values, which must be exactly 4.01, 7.00, or 10.01. The temperatures of these buffer solutions must be close to the same.

- ◆ If required, measure the temperature of the buffer solutions and ensure that the temperatures are close to the same during use of several buffer solutions.

**Material temperature differences between the buffer solutions reduce the quality of the calibration**

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**and therefore the subsequent measurement accuracy of the device.**

- ◆ Remove the protective cap from the sensor and rinse the sensor for approx. 10 minutes in a pH neutral liquid or in tap water.

### Calibration

- ◆ To carry out the calibration, briefly press the »CAL« key. This starts the calibration mode of the device.

**Calibration can be discontinued at any time. To do so, briefly press the »ON/HOLD« key.**

**Discontinuing calibration or an »ERROR« while determining the calibration values (see below) leads to lower measurement accuracy.**

**If the device does not receive any key inputs for several minutes while in calibration mode, it automatically switches back to measurement mode.**

- ◆ First, enter the temperature of the buffer solution. To do so, press the ▼ or ▲ keys until the actual temperature of the buffer solution(s) is shown on the display.
- ◆ Press the »CAL« key.

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


- ◆ Next, specify if you want to calibrate the device gradient using one, two, or three buffer solutions. To do so, press the ▼ or ▲, key until »1PC«, »2PC«, or »3PC« is shown.
- ◆ Press the »CAL« key. The display shows »CAL 1«.
- ◆ Insert the sensor into the first buffer solution and then press the »CAL« key again.

The device now determines the first calibration value. If the calibration value is valid, »OK« appears on the display after some time; otherwise »ERROR« appears (see note on page 26). In the latter case, repeat the process with another buffer solution.

If you previously selected »2PC« or »3PC«, »CAL 2« is now shown on the display. Otherwise, calibration is complete and the device switches back to measurement mode.

- ◆ If »CAL 2« is shown on the display, insert the sensor into the second buffer solution and press the »CAL« key.

**The order of the buffer solutions with the various pH values is random.**

The device now determines the second calibration value. If the calibration value is valid, »OK« appears on the display after some time;

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otherwise »ERROR« appears (see note on page 26). In the latter case, the device requests the second buffer solution (»CAL 2«) again.

**If applicable, use a buffer solution with a pH value that has not yet been used during this calibration session.**

If you previously selected »3PC«, »CAL 3« is now shown on the display. Otherwise, calibration is complete and the device switches back to measurement mode.

◆ If »CAL 3« is shown on the display, insert the sensor into the third buffer solution and press the »CAL« key.

The device now determines the third calibration value. If the calibration value is valid, »OK« appears on the display after some time; otherwise »ERROR« appears (see note on page 26). In the latter case, the device requests the third buffer solution (»CAL 3«) again.

**If applicable, use a buffer solution with a pH value that has not yet been used during this calibration session.**

If the device has determined a valid third calibration value, calibration is complete and the device switches back to measurement mode.

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
## What to do when...

the following is shown on the display:

Possible cause	Removal
Measurement range exceeded Sensor defective	Observe measurement range Call service
Below measurement range Short-circuit of the sensor	Observe measurement range Call service

### Changing the battery

The battery symbol in the display shows the battery condition.

Battery ok

Battery ok

Battery low

Battery change required

To change the lithium battery, the device must be opened. Proceed as follows:

**Use suitable precautions such as the use of a grounded wrist strap in order to prevent electrostatic discharge while changing the battery!**

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

***Electrostatic discharges can destroy the device!***

- ◆ First, carefully remove the two rubber stoppers (1) on the lower side of the device, e.g. using forceps.
- ◆ Now completely remove the two visible Phillips-head screws (2) with a suitable screwdriver (PZ1)
- ◆ Hold the bottom of the device (3) with one hand and carefully pull the upper part (4) off upwards with the other hand until the circuit board is exposed.

Now you can see the battery mount (7) with the battery (6) on the circuit board (5).

- ◆ Pull the spent battery from the battery mount in the direction of the arrow.
- ◆ Take the new battery with fat-free fingers and insert it into the battery mount. The plus sign on the battery must face upwards so it is visible.
- ◆ Reassemble the upper and lower parts of the housing.
- ◆ Tighten the two screws (only with moderate force, 0.4 Nm) and replace the rubber stoppers (the slightly beveled edges facing out).

**Dispose of the spent battery in an environmentally sound manner!**

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## Cleaning and care

Clean the device with a damp cloth.

**Do not use any solvents, e.g. acetone, since they attack synthetic materials.**

## Disposal

Should the device become unusable, you must dispose of it in a professional and environmentally sound manner.

**Do not dispose of the device in household garbage under any circumstances.**

**Dispose of the battery at a designated collection point.**

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## Technical data

Measurement range:	pH 0 to 14
Resolution of the measurement signal:	pH 0,01
Measurement accuracy of the electronics:	pH $\pm 0,03$
Temperature of the substance being measured:	
Adjustment range:	0 to 120°C (32 bis 248°F)
Increment:	0,5°C (0,5°F)
Operating temperature:	0 °C ... +50 °C (32 to 122°F)
Category temperature:	-10 °C ... +60 °C (14 to 140°F)
Storage temperature:	-25 °C ... +60 °C
Refresh rate:	1 to 15 seconds, adjustable
Automatic shutdown:	automatically after 2 hours, can be deactivated
Protection class:	IP 67
Dimensions: (L x W x H)	115 x 54 x 22 mm;
Housing material:	ABS
Weight	approx. 90 g
Sensor connector:	BNC
Battery	Lithium 3V/1Ah
Battery life:	up to 5 years

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Produktions- u. Vertriebs GmbH