

Data sheet

PRImeDETECT™ Legionella spp Plus Detection Kit (*Legionella* spp & *Legionella pneumophila*)

Cat. No: FP0025 (48 reactions)

Cat. No: FP0026 (96 reactions)

Introduction

Waterborne pathogens and related diseases are a major public health concern worldwide, not only by the morbidity and mortality that they cause, but by the high cost that represents their prevention and treatment. These diseases are directly related to environmental deterioration and pollution.

Legionella contamination in air conditioners and water supply systems poses a serious health concern. Canvax Biotech SL offers molecular microbiology PCR-based detection systems intended for the specific, rapid, and reliable detection of this pathogen in a user-friendly and cost-effective format.

PRImeDETECT™ Legionella spp Plus Detection Kit (*Legionella* spp & *Legionella pneumophila*) is based on amplification and detection of specific DNA fragments from *Legionella* spp and/or *Legionella pneumophila* by the real-time PCR method (multiplex PCR).

All reagents required for qPCR are provided ready to use as Multiplex PCR Master Mix (1 vial). The Multiplex PCR Master Mix contains the appropriate amounts of buffer, dNTPs, Hot-start DNA polymerase, DNA-free water and MgCl₂ to perform the number of reactions indicated in the kit. The Multiplex PCR Master Mix also includes an internal amplification control (IAC) whose detection indicates the absence of PCR inhibitors. Primers and probes for the amplification of IAC as well as for the amplification of the target gene are included in the Master Mix. The probe for the detection of target gene is labelled with the FAM (*Legionella* spp) and HEX fluorochrome (*L. pneumophila*), whereas the probe for the detection of IAC is labelled with the CY5 fluorochrome.

In addition, the kit includes positive control DNA and negative control. The positive control is supplied to demonstrate that the PCR amplification is working efficiently with the supplied components. To confirm absence of contamination, a negative control reaction should be included every time the kit is used

Include DNAREady lysis buffer to extract the DNA from the sample prior to PCR detection.

Each kit contains:

- ✓ Multiplex PCR Master Mix (1 vial)
- ✓ DNAREady lysis buffer (1 bottle)
- ✓ PCR Positive Control (1 vial)
- ✓ PCR Negative Control (1 vial)

Shipping and Storage

PRImeDETECT™ Legionella spp Plus Detection Kit is shipped at ambient temperature. On arrival the kit should be stored at -20°C. Avoid prolonged exposure to light. If stored correctly the kit will retain full activity for 12 months. The kit can be stored at 4°C for 1 month.

Technical features

- ✓ **Amplification limit:** *Legionella* spp: 5 UG per reaction. *L. pneumophila*: 100 UG per reaction.
- ✓ **Quantification limit:** *Legionella* spp: 50 UG/L. *L. pneumophila*: 1000 UG/L.
- ✓ **Quantification Dynamic range:** 5 logs
- ✓ **Simultaneously** detect and differentiate *Legionella* spp. and *L. pneumophila* in the same tube
- ✓ **Inclusivity:** Primers and probes have been tested satisfactorily in 14 strains of *L. pneumophila* and 30 species of the genus *Legionella*.
- ✓ **Compliance:** ISO/TS 12869:2012
- ✓ **Exclusivity:** No cross-reactions with other frequent bacterial agents have been observed.
- ✓ **Detection:** probe labelled with fluorescent dyes—*Legionella*: FAM-BHQ1a; *L. pneumophila*: HEX-BHQ1; IAC:CY5-BHQ2.
- ✓ **Thermal cycler:** Agilent Mx3005P, Applied Biosystems 7300, 7500 and other cyclers.

Distributed by:

Lab Unlimited
CARL STUART GROUP

Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432

Fax: (01) 4523967

E-mail: info@labunlimited.com

Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30

Fax: 08452 30 50 30

E-mail: info@labunlimited.co.uk

Web: www.labunlimited.co.uk

(continued on reverse side)

 **canvax**

PCR cycling conditions

Step	Time	Temperature
Initial denaturation	10 min	95°C
45 Cycles:	15 sec	95°C
	1 min	60°C
Melt analysis	Refer to instrument instructions	

Analysis of results

Follow instrument software instructions to generate cycle threshold (Ct) values from the acquired data. The user may also, optionally, analyze the melt profile of each reaction.

The quantity of DNA target in each sample can be calculated by referring to the positive control template Ct value.

PRODUCT USE LIMITATION

Not For Medical Diagnostic Use.

Distributed by:

Lab Unlimited
CARL STUART GROUP

Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432
Fax: (01) 4523967
E-mail: info@labunlimited.com
Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30
Fax: 08452 30 50 30
E-mail: info@labunlimited.co.uk
Web: www.labunlimited.co.uk