



titan

continuous processing
scale-up

Distributed by:

Lab Unlimited
CARL STUART GROUP

Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432
Fax: (01) 4523967
E-mail: info@labunlimited.com
Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30
Fax: 08452 30 50 30
E-mail: info@labunlimited.co.uk
Web: www.labunlimited.co.uk

» to manufacturing



Manufacturing

Using the same continuous processing modules available at lab and pilot scale, you can move straight into manufacturing scale with minimal fuss.

Titan systems can produce up to a tonne (process dependent) of product per day.

Titan systems are designed for long term operation with built-in safety features.

» specifications

- Reactor volumes:** 64 mL and 250 mL reactors can be added in series to cover a wide volume range
- Production rates:** Output suitable for laboratory, pilot plant and manufacturing environments
- System flow rates:** 1 mL/min to 1 L/min
- Pump flow rates:** 1 mL/min to 250 mL/min
- Pressures:** Up to 20 bar
- Temperature range:** -40°C to 250°C
- Types of processes:** Suited for both fast reactions with high mass transfer and heat transfer requirements, and slow reactions. Designed for use with multiple phases
- Chemical compatibility:** Glass, fluoropolymers, Hastelloy, and Elgiloy wetted parts
- Automation:** Control process parameters from the modules, via Titan software, or via pre-existing PLC or plant software
- Safety:** Interlocks, anti-tamper, pressure detection, and auto shut down features