


**A 101**

Upper basket  
 Height adjustable, for the optimum loading of inserts.



- Vertical clearance 160 +/- 30 mm, width 475 mm, depth 443 mm
- With integrated spray arm



## Product details - A 101

EAN: 4002515460627 / Article no. 69610100D / Mat. no. 09862340

Product category <span style="float: right;">^</span>	
Fresh-water dishwashers	•
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine/accessory categorisation <span style="float: right;">^</span>	
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8592	•
PG 8593	•
Product category <span style="float: right;">^</span>	
Carrier/basket carrier	•
Upper basket	•
Product properties <span style="float: right;">^</span>	
Material	Stainless steel
Colour	Stainless steel
Area of application <span style="float: right;">^</span>	
Reprocessing laboratory glassware	•
Reprocessing GYN instruments	•

Reprocessing ENT instruments	•
Reprocessing basins	•
Reprocessing kidney dishes	•
Reprocessing OP instruments	•
Reprocessing of different utensils	•

Dimensions and weight		^
External dimensions, height in mm	206	
External dimensions, width in mm	528	
External dimensions, depth in mm	527	
Net weight in kg	3,000	

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:

**Lab Unlimited**  
CARL STUART GROUP

Tallaght Business Park  
Whitestown, Dublin 24,  
Ireland  
D24 RFK3

Tel: (01) 4523432  
Fax: (01) 4523967  
E-mail: [info@labunlimited.com](mailto:info@labunlimited.com)  
Web: [www.labunlimited.com](http://www.labunlimited.com)

Quatro House, Frimley Road,  
Camberley,  
United Kingdom  
GU16 7ER

Tel: 08452 30 40 30  
Fax: 08452 30 50 30  
E-mail: [info@labunlimited.co.uk](mailto:info@labunlimited.co.uk)  
Web: [www.labunlimited.co.uk](http://www.labunlimited.co.uk)