



A 12/1

Insert

For the gentle cleaning of delicate items by interrupting the spray jet.

- With perforations 7 x 7 mm
- For use in upper and lower basket



Product details - A 12/1

EAN: 4002513975697 / Article no. 69501203D / Mat. no. 06629670

Product category		^
Washer-disinfectors, laboratory		•
Machine/accessory categorisation		^
PG 8535		•
PG 8536		•
PG 8583		•
PG 8583 CD		•
PG 8593		•
Product category		^
Retaining frame		•
Product properties		^
Material	Stainless steel	
Colour	White	
Area of application		^
Reprocessing laboratory glassware		•
Dimensions and weight		^
External dimensions, width in mm	224	
External dimensions, depth in mm	429	
Net weight in kg	0,372	

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:



Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432
Fax: (01) 4523967
E-mail: info@labunlimited.com
Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30
Fax: 08452 30 50 30
E-mail: info@labunlimited.co.uk
Web: www.labunlimited.co.uk