



A 13

Lid







For fixing delicate wash items.

- Lid in stainless steel
- 1 mm wire, 8 mm mesh size, 4 mm frames



Product details - A 13

EAN: 4002511117204 / Article no. 69745001D / Mat. no. 03810200

Product category 	
Washer-disinfectors, laboratory	•
Machine/accessory categorisation 	
PG 8535	•
PG 8536	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category 	
Lid	•
Product properties 	
Material	Stainless steel
Colour	Stainless steel
Area of application 	
Reprocessing laboratory glassware	•
Dimensions and weight 	
External dimensions, height in mm	20
External dimensions, width in mm	200
External dimensions, depth in mm	320
Net weight in kg	0,135

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:



Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432
Fax: (01) 4523967
E-mail: info@labunlimited.com
Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30
Fax: 08452 30 50 30
E-mail: info@labunlimited.co.uk
Web: www.labunlimited.co.uk