

### A 151

Lower basket  
For holding standard DIN mesh trays as well as various inserts.



- Vertical clearance dependent on upper basket selected ⓘ
- Loading width: 490 mm, depth: 435 mm



## Product details - A 151

EAN: 4002515460689 / Article no. 69615100D / Mat. no. 09862410

|                                  |                 |   |
|----------------------------------|-----------------|---|
| Product category                 |                 | ^ |
| Washer-disinfectors, laboratory  | •               |   |
| Thermal disinfectors             | •               |   |
| Washer-disinfectors, medical     | •               |   |
| Machine/accessory categorisation |                 | ^ |
| PG 8581                          | •               |   |
| PG 8562                          | •               |   |
| PG 8582                          | •               |   |
| PG 8582 CD                       | •               |   |
| PG 8504                          | •               |   |
| PG 8583                          | •               |   |
| PG 8583 CD                       | •               |   |
| PG 8591                          | •               |   |
| PG 8592                          | •               |   |
| PG 8593                          | •               |   |
| Product category                 |                 | ^ |
| Carrier/basket carrier           | •               |   |
| Lower basket                     | •               |   |
| Product properties               |                 | ^ |
| Material                         | Stainless steel |   |
| Colour                           | Stainless steel |   |
| Area of application              |                 | ^ |

|                                    |   |
|------------------------------------|---|
| Reprocessing laboratory glassware  | • |
| Reprocessing GYN instruments       | • |
| Reprocessing ENT instruments       | • |
| Reprocessing dental instruments    | • |
| Reprocessing theatre shoes         | • |
| Reprocessing basins                | • |
| Reprocessing kidney dishes         | • |
| Reprocessing baby bottles          | • |
| Reprocessing OP instruments        | • |
| Reprocessing of mesh trays         | • |
| Reprocessing of different utensils | • |

|          |  |   |
|----------|--|---|
| Capacity |  | ^ |
|----------|--|---|

|                          |   |
|--------------------------|---|
| DIN mesh trays [number]  | 2 |
| SPRI mesh trays [number] | 1 |

|                       |  |   |
|-----------------------|--|---|
| Dimensions and weight |  | ^ |
|-----------------------|--|---|

|                                   |       |
|-----------------------------------|-------|
| External dimensions, height in mm | 88    |
| External dimensions, width in mm  | 529   |
| External dimensions, depth in mm  | 522   |
| Net weight in kg                  | 1,800 |

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:

**Lab Unlimited**  
CARL STUART GROUP

Tallaght Business Park  
Whitestown, Dublin 24,  
Ireland  
D24 RFK3

Tel: (01) 4523432  
Fax: (01) 4523967

E-mail: [info@labunlimited.com](mailto:info@labunlimited.com)  
Web: [www.labunlimited.com](http://www.labunlimited.com)

Quatro House, Frimley Road,  
Camberley,  
United Kingdom  
GU16 7ER

Tel: 08452 30 40 30  
Fax: 08452 30 50 30

E-mail: [info@labunlimited.co.uk](mailto:info@labunlimited.co.uk)  
Web: [www.labunlimited.co.uk](http://www.labunlimited.co.uk)