


A 304

Module

For the optimum loading of up to 98 vials.

- Suitable for 98 centrifuge tubes, phials, test tubes, for example
- The loading heights differ depending on the positioning ⓘ



Product details - A 304

EAN: 4002515460832 / Article no. 69630400D / Mat. no. 09862570

Product category		^
Washer-disinfectors, laboratory		•
Machine/accessory categorisation		^
PG 8504		•
PG 8583		•
PG 8583 CD		•
PG 8593		•
Product category		^
Insert/module for test tubes		•
Product properties		^
Material	Stainless steel Plastic	
Colour	Stainless steel Black	
Area of application		^
Reprocessing pipettes/vials		•
Capacity		^
Vials [number]	98	
Dimensions and weight		^
External dimensions, height in mm	130	
External dimensions, width in mm	222	
External dimensions, depth in mm	471	
Net weight in kg	2,300	

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:



Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432
Fax: (01) 4523967
E-mail: info@labunlimited.com
Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30
Fax: 08452 30 50 30
E-mail: info@labunlimited.co.uk
Web: www.labunlimited.co.uk