

Accessories for all Chambers and Glove Boxes: Glove Box & Isolator Glove Specifications

We offer 22 variations of gauntlet gloves and replacement hands for glove boxes and isolators. All of the gauntlet gloves (sleeves) can also be mounted onto our conversion kits. Conversion kits (formed plastic caps) can be mounted onto any of our glove boxes or isolators.

Specification Chart

| Product Number | Description | Material | Port Diameter | Hand Size | Thickness |
|--|-----------------------------|----------|---------------|-----------|-----------|
| 800-GH/S | Sleeve Only | Hypalon | 8" | n/a | .015" |
| 800-GOS | Glove over sleeve connector | | | | |
| Optional hand replacement for use with 800-GH/S, sleeve only. 800-GOS, glove over sleeve connector, is needed for use with the hand replacements and sleeves. | | | | | |
| 800-GA/BUT/7 | Hand Replacement | Butyl | n/a | 7" | .015" |
| 800-GA/BUT/8 | Hand Replacement | Butyl | n/a | 8" | .015" |
| 800-GA/BUT/9 | Hand Replacement | Butyl | n/a | 9" | .015" |
| 800-GA/BUT/10 | Hand Replacement | Butyl | n/a | 10" | .015" |
| 800-GA/HYP/7 | Hand Replacement | Hypalon | n/a | 7" | .015" |
| 800-GA/HYP/8 | Hand Replacement | Hypalon | n/a | 8" | .015" |
| 800-GA/HYP/9 | Hand Replacement | Hypalon | n/a | 9" | .015" |
| 800-GA/HYP/10 | Hand Replacement | Hypalon | n/a | 10" | .015" |
| 800-GA/NEO/7 | Hand Replacement | Neoprene | n/a | 7" | .015" |
| 800-GA/NEO/8 | Hand Replacement | Neoprene | n/a | 8" | .015" |
| 800-GA/NEO/9 | Hand Replacement | Neoprene | n/a | 9" | .015" |
| 800-GA/NEO/10 | Hand Replacement | Neoprene | n/a | 10" | .015" |
| 800-GA/NIT/7 | Hand Replacement | Nitrile | n/a | 7" | .015" |
| 800-GA/NIT/8 | Hand Replacement | Nitrile | n/a | 8" | .015" |
| 800-GA/NIT/9 | Hand Replacement | Nitrile | n/a | 9" | .015" |
| 800-GA/NIT/10 | Hand Replacement | Nitrile | n/a | 10" | .015" |
| One-Piece Gloves | | | | | |
| 800-GB | One-Piece Glove | Butyl | 8" | 9" | .015" |
| 800-GB/7 | One-Piece Glove | Butyl | 8" | 7" | .015" |
| 800-GB/8 | One-Piece Glove | Butyl | 8" | 8" | .015" |
| 800-GB/9/6 | One-Piece Glove | Butyl | 7" | 9" | .015" |
| 800-GB/9/10 | One-Piece Glove | Butyl | 10" | 9" | .015" |
| 800-GH | One-Piece Glove | Hypalon | 8" | 9" | .015" |
| 800-GH/7 | Two-Piece Glove | Hypalon | 8" | 7" | .015" |
| 800-GH/8 | One-Piece Glove | Hypalon | 8" | 8" | .015" |
| 800-GH/10 | One-Piece Glove | Hypalon | 8" | 10" | .015" |
| 800-GH/8/10 | One-Piece Glove | Hypalon | 10" | 8" | .015" |
| 800-GH/9/10 | One-Piece Glove | Hypalon | 10" | 9" | .015" |
| 800-GH/9/6 | One-Piece Glove | Hypalon | 7" | 9" | .015" |
| 815-PGB, "La Petite" Glove Box Gloves | | | | | |
| 800-GLLP | One-Piece Glove | Latex | 7" | 9" | .015" |
| 800-GLLP/7 | One-Piece Glove | Latex | 7" | 7" | .015" |
| 800-GLLP/10 | One-Piece Glove | Latex | 7" | 10" | .015" |
| 800-GN | One-Piece Glove | Neoprene | 8" | 9" | .015" |
| 800-GN/8 | One-Piece Glove | Neoprene | 8" | 8" | .015" |
| 800-GN/10 | One-Piece Glove | Neoprene | 8" | 10" | .015" |
| 800-GN/11 | One-Piece Glove | Neoprene | 8" | 11" | .015" |
| 800-GN/9/6 | One-Piece Glove | Neoprene | 7" | 9" | .015" |