



Pipettor Compatibility

Brand	Model	Volume Range	Channels	6.283 565/7 0.1-10µl	6.283 568/9 200µl	6.283 570/2 1250µl
LLG Labware	proMLP	0.2-2 µl	1	x		
LLG Labware	proMLP	0.5-10 µl	1	x		
LLG Labware	proMLP	2-20 µl	1		x	
LLG Labware	proMLP	10-100 µl	1		x	
LLG Labware	proMLP	20-200 µl	1		x	
LLG Labware	proMLP	100-1000 µl	1			x
LLG Labware	Microliter pipette, fix	5 µl	1		x	
LLG Labware	Microliter pipette, fix	10 µl	1		x	
LLG Labware	Microliter pipette, fix	20 µl	1		x	
LLG Labware	Microliter pipette, fix	25 µl	1		x	
LLG Labware	Microliter pipette, fix	50 µl	1		x	
LLG Labware	Microliter pipette, fix	100 µl	1		x	
LLG Labware	Microliter pipette, fix	200 µl	1		x	
LLG Labware	Microliter pipette, fix	250 µl	1			x
LLG Labware	Microliter pipette, fix	500 µl	1			x
LLG Labware	Microliter pipette, fix	1000 µl	1			x
LLG Labware	Microliter pipette, variable	0.1-2.5 µl	1	x		
LLG Labware	Microliter pipette, variable	0.5-10 µl	1	x		
LLG Labware	Microliter pipette, variable	2-20 µl	1	x		
LLG Labware	Microliter pipette, variable	10-100 µl	1		x	
LLG Labware	Microliter pipette, variable	20-200 µl	1		x	
LLG Labware	Microliter pipette, variable	100-1000 µl	1			x
Biohit	mLine	0.1-3 µl, 0.5-10 µl	1	x		
Biohit	mLine	2-20 µl, 10-100 µl, 20-200 µl	1		x	
Biohit	mLine	100-1000 µl	1			x
Biohit	mLine	0.5-10 µl	8/12	x		
Biohit	mLine	5-100 µl	8/12		x	
Biohit	mLine	30-300 µl	8/12			
Biohit	Proline	0.1-2.5 µl, 0.5-10 µl	1	x		
Biohit	Proline	2-20 µl, 5-50 µl, 10-100 µl, 20-200 µl	1		x	
Biohit	Proline	100-1000 µl	1			x
Biohit	Proline	0.5-10 µl	8/12	x		
Biohit	Proline	5-50 µl	8/12		x	
Biohit	Proline	50-300 µl	8/12			
Brand	Transferpette	0.1-1 µl, 0.5-10 µl	1	x		
Brand	Transferpette	2-20 µl micro	1	x		
Brand	Transferpette	2-20 µl, 5-50 µl, 10-100 µl, 20-200 µl	1		x	
Brand	Transferpette	100-1000 µl	1			x
Brand	Transferpette	0.5-10 µl	8/12	x		
Brand	Transferpette	2-20 µl	8/12	x		
Brand	Transferpette	2.5-25 µl, 5-50 µl, 10-100 µl, 20-200 µl	8/12		x	
Brand	Transferpette	30-300 µl	8/12			
Brand	Transferpette S	0.1-1 µl, 0.1-2.5 µl, 0.5-10 µl	1	x		
Brand	Transferpette S	2-20 µl, 10-100 µl, 20-200 µl	1		x	
Brand	Transferpette S	100-1000 µl	1			x
Brand	Transferpette S	0.5-10 µl	8/12	x		
Brand	Transferpette S	5-50 µl, 10-100 µl, 20-200 µl	8/12		x	
Brand	Transferpette S	30-300 µl	8/12			
Eppendorf	Research Plus	0.1-2.5 ul, 0.5-10 µl	1	x		
Eppendorf	Research Plus	2-20 µl micro	1	x		



LLG Pipette Tips **ULTRALOW**

Pipettor Compatibility

Brand	Model	Volume Range	Channels	6.283 565/7 0.1-10µl	6.283 568/9 200µl	6.283 570/2 1250µl
Eppendorf	Research Plus	2-20 µl, 10-100 µl, 20-200 µl	1		x	
Eppendorf	Research Plus	30-300 µl	1			
Eppendorf	Research Plus	100-1000 µl	1			x
Eppendorf	Research Plus	0.5-10 µl	8/12	x		
Eppendorf	Research Plus	10-100 µl	8/12		x	
Eppendorf	Research Plus	30-300 µl	8/12			
Eppendorf	Reference	0.1-2.5 ul, 0.5-10 µl	1	x		
Eppendorf	Reference	2-20 µl micro	1	x		
Eppendorf	Reference	2-20 µl, 10-100 µl, 50-200 µl	1		x	
Eppendorf	Reference	100-1000 µl	1			x
Gilson	Pipetman Classic, G, L: P2	0.2-2 ul, 1-10 ul	1	x		
Gilson	Pipetman Classic, G, L:	2-20 ul, 20-100 ul, 50-200 ul	1		x	
Gilson	Pipetman Classic, G, L:	200-1000 ul	1			x
Gilson	Pipetman P10L	0.5-10 µl, 2-20 µl, 20-200 µl	8/12		x	
Gilson	Pipetman P300L	20-300 µl	8/12			
Nichiryo	Nichipet Premium, EX, EX	0.1-2 µl	1	x		
Nichiryo	Nichipet Premium, EX, EX	2-20 µl, 10-100 µl, 20-200 µl	1		x	
Nichiryo	Nichipet Premium, EX, EX	100-1000 µl	1			x
Nichiryo	Nichipet 7000	0.5-10 µl	8/12	x		
Nichiryo	Nichipet 7000	5-50 µl, 40-200 µl	8/12		x	
Nichiryo	Nichipet 7000	50-300 µl	8/12			
Rainin	Pipet-Lite SL-2XLS	0.1-2 µl, 0.5-10 µl	1	x		
Rainin	Pipet-Lite SL-20XLS	2-20 µl, 10-100 µl, 20-200 µl	1		x	
Rainin	Pipet-Lite	20-300 µl	1			
Rainin	Pipet-Lite	100-1000 µl	1			x
Rainin	Classic PR-2	0.1-2 µl, 0.5-10 µl	1	x		
Rainin	Classic PR-20	2-20 µl, 10-100 µl, 20-200 µl	1		x	
Rainin	Classic PR-1000	100-1000 µl	1			x
Socorex	Calibra® 822	0.2-2 µl, 1-10 µl	1	x		
Socorex	Calibra® 822	1-10 µl, 2-20 µl, 10-100 µl, 20-200 µl	1		x	
Socorex	Calibra® 822	100-1000 µl	1			x
Socorex	Calibra® 852	1-10 µl	8	x		
Socorex	Calibra® 852	10-100 µl, 20-200 µl	8/12		x	
Socorex	Acura 825	0.1-2 µl, 0.5-10 µl	1	x		
Socorex	Acura 825	1-10 µl, 2-20 µl, 5-50 µl 10-100 µl, 20-200 µl	1		x	
Socorex	Acura 825	100-1000 µl	1			x
Socorex	Acura 855	0.5-10 µl	8/12	x		
Socorex	Acura 855	5-50 µl, 2-20 µl	8/12		x	
Socorex	Acura 855	40-350 µl	8/12			
Thermo/Finn	Finnpette 4500	0.2-2 µl, micro, 0.5-5 µl	1	x		
Thermo/Finn	Finnpette 4500	0.5-5 µl, 2-20 µl, 5-50 µl	1		x	
Thermo/Finn	Finnpette 4500	100-1000 µl	1			x
Thermo/Finn	Finnpette 4500	0.5-10 µl	8/12	x		
Thermo/Finn	Finnpette 4500	5-50 µl	8/12		x	
Thermo/Finn	Finnpette 4500	50-300 µl	8/12			