


## E 111

## Insert

For the optimum loading of beakers from 600 to 1,000 ml.

- Suitable for 8 glass beakers between 600 - 1000 ml
- equipped with 8 x 3 stay bars
- For use in lower basket



## Product details - E 111

EAN: 4002511115743 / Article no. 69511101D / Mat. no. 03808420

Product category <span style="float: right;">^</span>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine/accessory categorisation <span style="float: right;">^</span>	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category <span style="float: right;">^</span>	
Insert/module for beakers	•
Product properties <span style="float: right;">^</span>	
Material	Stainless steel
Colour	Stainless steel
Area of application <span style="float: right;">^</span>	
Reprocessing laboratory glassware	•
Dimensions and weight <span style="float: right;">^</span>	
External dimensions, height in mm	205
External dimensions, width in mm	230
External dimensions, depth in mm	460

Net weight in kg

0,960

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:

**Lab Unlimited**  
CARL STUART GROUP

Tallaght Business Park  
Whitestown, Dublin 24,  
Ireland  
D24 RFK3

Tel: (01) 4523432

Fax: (01) 4523967

E-mail: [info@labunlimited.com](mailto:info@labunlimited.com)

Web: [www.labunlimited.com](http://www.labunlimited.com)

Quatro House, Frimley Road,  
Camberley,  
United Kingdom  
GU16 7ER

Tel: 08452 30 40 30

Fax: 08452 30 50 30

E-mail: [info@labunlimited.co.uk](mailto:info@labunlimited.co.uk)

Web: [www.labunlimited.co.uk](http://www.labunlimited.co.uk)