

**E 137**

Insert  
for optimum loading of 56 half petri dishes.

- Suitable for 56 half petri dishes with Ø 100 mm
- equipped with 56 holders, height 70 mm, distance apart approx. 26 mm
- For use in lower basket



## Product details - E 137

EAN: 4002511115958 / Article no. 69513701D / Mat. no. 03830290

Product category 	
Washer-disinfectors, laboratory	•
Machine/accessory categorisation 	
PG 8535	•
PG 8536	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category 	
Insert/module for petri dishes	•
Product properties 	
Material	Stainless steel
Colour	Stainless steel
Distance between rails in mm	26
Rail height in mm	70
Area of application 	
Reprocessing laboratory glassware	•
Capacity 	
Petri dishes [number]	56
Dimensions and weight 	
External dimensions, height in mm	485

External dimensions, width in mm	445
External dimensions, depth in mm	95
Net weight in kg	2,100

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:

**Lab Unlimited**  
CARL STUART GROUP

Tallaght Business Park  
Whitestown, Dublin 24,  
Ireland  
D24 RFK3

Tel: (01) 4523432  
Fax: (01) 4523967

E-mail: [info@labunlimited.com](mailto:info@labunlimited.com)  
Web: [www.labunlimited.com](http://www.labunlimited.com)

Quatro House, Frimley Road,  
Camberley,  
United Kingdom  
GU16 7ER

Tel: 08452 30 40 30  
Fax: 08452 30 50 30

E-mail: [info@labunlimited.co.uk](mailto:info@labunlimited.co.uk)  
Web: [www.labunlimited.co.uk](http://www.labunlimited.co.uk)