



E 139/1

Insert

For the optimum loading of test tubes up to 12 x 200 mm.

- Suitable for test tubes up to 12 x 200 mm
- Divided into 6 surfaces
- Mesh width 9 x 9 mm
- For use in lower basket
- Including lid A 13



Product details - E 139/1

EAN: 4002514231082 / Article no. 69513902D / Mat. no. 06907660

Product category ^	
Washer-disinfectors, laboratory	•
Machine/accessory categorisation ^	
PG 8535	•
PG 8536	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category ^	
Insert/module for test tubes	•
Product properties ^	
Material	Stainless steel
Colour	Stainless steel
Number of compartments	6
Area of application ^	
Reprocessing laboratory glassware	•
Capacity ^	
Test tubes [number]	200
Dimensions and weight ^	
External dimensions, height in mm	243
External dimensions, width in mm	200

External dimensions, depth in mm

320

Net weight in kg

1,260

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:

Lab Unlimited
CARL STUART GROUP

Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432

Fax: (01) 4523967

E-mail: info@labunlimited.com

Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30

Fax: 08452 30 50 30

E-mail: info@labunlimited.co.uk

Web: www.labunlimited.co.uk