



E 329

Mobile unit
for the optimum loading of narrow necked
glassware.

- equipped with 39 injector nozzles



Product details - E 329

EAN: 4002511116108 / Article no. 69532901D / Mat. no. 03809340

Product category		^
Washer-disinfectors, laboratory		•
Machine/accessory categorisation		^
PG 8535		•
PG 8536		•
Product category		^
Mobile unit for narrow-necked glassware		•
Product properties		^
Material	Stainless steel	
Colour	Stainless steel	
Area of application		^
Reprocessing laboratory glassware		•
Capacity		^
Narrow-necked glassware [number]	39	
Dimensions and weight		^
External dimensions, height in mm	502	
External dimensions, width in mm	535	
External dimensions, depth in mm	517	
Net weight in kg	5,500	
Standard accessories		^
ID 90 injector nozzles Ø 2.5 mm length 90 mm [number]	4	

ID 110, injector nozzles Ø 2.5 mm, length 110 mm [number]	5
ID 140, injector nozzles Ø 4 mm, length 140 mm [number]	5
ID 160, injector nozzles Ø 4 mm, length 160 mm [number]	3 5
ID 180, injector nozzles Ø 4 mm, length 180 mm [number]	5
ID 200, injector nozzles Ø 6 mm, length 200 mm [number]	5
ID 220, injector nozzles Ø 6 mm, length 220 mm [number]	5
ID 240, injector nozzles Ø 6 mm, length 240 mm [number]	5
Rinse nozzle for powder dispenser [number]	1

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:

Lab Unlimited
CARL STUART GROUP

Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432

Fax: (01) 4523967

E-mail: info@labunlimited.com

Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30

Fax: 08452 30 50 30

E-mail: info@labunlimited.co.uk

Web: www.labunlimited.co.uk