



E 353

Clamps

Height-adjustable for optimum positioning of laboratory glassware on the nozzle.

- Suitable for injector nozzle E 351
- Height adjustable locking
- Suitable for bottles, Erlenmeyer flasks and round bottomed flasks for example



Product details - E 353

EAN: 4002511117099 / Article no. 69735301D / Mat. no. 03809530

Product category		^
Washer-disinfectors, laboratory		•
Large-capacity washer-disinfectors, laboratory		•
Machine/accessory categorisation		^
G 7825		•
PG 8527		•
PG 8535		•
PG 8536		•
PG 8583		•
PG 8583 CD		•
PG 8593		•
Product category		^
Other accessories		•
Product properties		^
Material	Stainless steel Wire, rilsan-coated	
Colour	Stainless steel White	
Area of application		^
Reprocessing laboratory glassware		•
Dimensions and weight		^
External dimensions, depth in mm		160
Net weight in kg		0,040

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:



Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432
Fax: (01) 4523967
E-mail: info@labunlimited.com
Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30
Fax: 08452 30 50 30
E-mail: info@labunlimited.co.uk
Web: www.labunlimited.co.uk