



## E 380

Mobile unit  
for the optimum loading of narrow necked  
glassware.

- equipped with 32 nozzles and locking mechanisms
- Equipped with connection for hot air drying
- Bracket for ML/2 for automatic mobile unit recognition



## Product details - E 380

EAN: 4002511116467 / Article no. 69538001D / Mat. no. 03809880

Product category <span style="float: right;">^</span>	
Washer-disinfectors, laboratory	•
Machine/accessory categorisation <span style="float: right;">^</span>	
PG 8535	•
PG 8536	•
Product category <span style="float: right;">^</span>	
Mobile unit with connection for hot air drying	•
Mobile unit for narrow-necked glassware	•
Product properties <span style="float: right;">^</span>	
Material	Stainless steel
Colour	Stainless steel
Area of application <span style="float: right;">^</span>	
Reprocessing laboratory glassware	•
Capacity <span style="float: right;">^</span>	
Narrow-necked glassware [number]	32
Dimensions and weight <span style="float: right;">^</span>	
External dimensions, height in mm	502
External dimensions, width in mm	535
External dimensions, depth in mm	517
Net weight in kg	7,900
Standard accessories <span style="float: right;">^</span>	

Standard accessories	
E 351 injector nozzles Ø 4 mm, length 160 mm [number]	14
E 352 injector nozzles Ø 6 mm, length 220 mm [number]	18
E 353 clips for injector nozzles E 351 [number]	14
E 354 clips for injector nozzles E 352 [number]	18
Rinse nozzle for powder dispenser [number]	1

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:

**Lab Unlimited**  
CARL STUART GROUP

Tallaght Business Park  
Whitestown, Dublin 24,  
Ireland  
D24 RFK3

Tel: (01) 4523432

Fax: (01) 4523967

E-mail: [info@labunlimited.com](mailto:info@labunlimited.com)

Web: [www.labunlimited.com](http://www.labunlimited.com)

Quatro House, Frimley Road,  
Camberley,  
United Kingdom  
GU16 7ER

Tel: 08452 30 40 30

Fax: 08452 30 50 30

E-mail: [info@labunlimited.co.uk](mailto:info@labunlimited.co.uk)

Web: [www.labunlimited.co.uk](http://www.labunlimited.co.uk)