



E 402

Insert
for optimum loading of up to 44 watch glasses.

- Suitable for 44 watch glasses Ø 80 - 125 mm
- equipped with 23 bars spaced 15 mm apart



Product details - E 402

EAN: 4002511116627 / Article no. 69540201D / Mat. no. 03830420

Product category ^	
Washer-disinfectors, laboratory	•
Machine/accessory categorisation ^	
PG 8535	•
PG 8536	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category ^	
Insert/module for watch glasses	•
Product properties ^	
Material	Stainless steel
Colour	Stainless steel
Distance between rails in mm	15
Area of application ^	
Reprocessing laboratory glassware	•
Capacity ^	
Watch glasses [number]	44
Dimensions and weight ^	
External dimensions, height in mm	53
External dimensions, width in mm	200

External dimensions, depth in mm

445

Net weight in kg

0,910

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:

Lab Unlimited
CARL STUART GROUP

Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432
Fax: (01) 4523967

E-mail: info@labunlimited.com
Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30
Fax: 08452 30 50 30

E-mail: info@labunlimited.co.uk
Web: www.labunlimited.co.uk