

## E 403

Insert  
for optimum loading of up to 105 watch glasses.

- Suitable for 105 watch glasses Ø 50 - 60 mm
- equipped with 36 bars spaced 9 mm apart



## Product details - E 403

EAN: 4002511116634 / Article no. 69540301D / Mat. no. 03830430

Product category		^
Washer-disinfectors, laboratory		•
Product category		^
Insert/module for watch glasses		•
Product properties		^
Material	Stainless steel	
Colour	Stainless steel	
Distance between rails in mm	9	
Area of application		^
Reprocessing laboratory glassware		•
Capacity		^
Watch glasses [number]	105	
Dimensions and weight		^
External dimensions, height in mm	35	
External dimensions, width in mm	445	
External dimensions, depth in mm	200	
Net weight in kg	1,450	

Distributed by:



Tallaght Business Park  
Whitestown, Dublin 24,  
Ireland  
D24 RFK3

Tel: (01) 4523432  
Fax: (01) 4523967  
E-mail: [info@labunlimited.com](mailto:info@labunlimited.com)  
Web: [www.labunlimited.com](http://www.labunlimited.com)

Quatro House, Frimley Road,  
Camberley,  
United Kingdom  
GU16 7ER

Tel: 08452 30 40 30  
Fax: 08452 30 50 30  
E-mail: [info@labunlimited.co.uk](mailto:info@labunlimited.co.uk)  
Web: [www.labunlimited.co.uk](http://www.labunlimited.co.uk)