



E 406

Mobile unit
for the optimum loading of up to 116 pipettes.

- Suitable for 116 pipettes up to 450 mm
- Surface area 16 x 16 mm
- Complete with frame height 150 mm
- Without bracket for magnetic strip for mobile unit recognition



Product details - E 406

EAN: 4002511116665 / Article no. 69540601D / Mat. no. 03810110

Product category		^
Washer-disinfectors, laboratory		•
Machine/accessory categorisation		^
PG 8535		•
PG 8536		•
Product category		^
Mobile unit for pipettes		•
Product properties		^
Material	Stainless steel Plastic	
Colour	Stainless steel Black	
Area of application		^
Reprocessing pipettes/vials		•
Capacity		^
Pipettes [number]	116	
Dimensions and weight		^
External dimensions, height in mm	502	
External dimensions, width in mm	535	
External dimensions, depth in mm	517	
Net weight in kg	4,900	

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:



Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432
Fax: (01) 4523967
E-mail: info@labunlimited.com
Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30
Fax: 08452 30 50 30
E-mail: info@labunlimited.co.uk
Web: www.labunlimited.co.uk