

# Accessories for all Chambers and Glove Boxes: Glove Box & Isolator Glove Specifications

We offer 22 variations of gauntlet gloves and replacement hands for glove boxes and isolators. All of the gauntlet gloves (sleeves) can also be mounted onto our conversion kits. Conversion kits (formed plastic caps) can be mounted onto any of our glove boxes or isolators.

## Specification Chart

Product Number	Description	Material	Port Diameter	Hand Size	Thickness
800-GH/S	Sleeve Only	Hypalon	8"	n/a	.015"
800-GOS	Glove over sleeve connector				
<b>Optional hand replacement for use with 800-GH/S, sleeve only. 800-GOS, glove over sleeve connector, is needed for use with the hand replacements and sleeves.</b>					
800-GA/BUT/7	Hand Replacement	Butyl	n/a	7"	.015"
800-GA/BUT/8	Hand Replacement	Butyl	n/a	8"	.015"
800-GA/BUT/9	Hand Replacement	Butyl	n/a	9"	.015"
800-GA/BUT/10	Hand Replacement	Butyl	n/a	10"	.015"
800-GA/HYP/7	Hand Replacement	Hypalon	n/a	7"	.015"
800-GA/HYP/8	Hand Replacement	Hypalon	n/a	8"	.015"
800-GA/HYP/9	Hand Replacement	Hypalon	n/a	9"	.015"
800-GA/HYP/10	Hand Replacement	Hypalon	n/a	10"	.015"
800-GA/NEO/7	Hand Replacement	Neoprene	n/a	7"	.015"
800-GA/NEO/8	Hand Replacement	Neoprene	n/a	8"	.015"
800-GA/NEO/9	Hand Replacement	Neoprene	n/a	9"	.015"
800-GA/NEO/10	Hand Replacement	Neoprene	n/a	10"	.015"
800-GA/NIT/7	Hand Replacement	Nitrile	n/a	7"	.015"
800-GA/NIT/8	Hand Replacement	Nitrile	n/a	8"	.015"
800-GA/NIT/9	Hand Replacement	Nitrile	n/a	9"	.015"
800-GA/NIT/10	Hand Replacement	Nitrile	n/a	10"	.015"
<b>One-Piece Gloves</b>					
800-GB	One-Piece Glove	Butyl	8"	9"	.015"
800-GB/7	One-Piece Glove	Butyl	8"	7"	.015"
800-GB/8	One-Piece Glove	Butyl	8"	8"	.015"
800-GB/9/6	One-Piece Glove	Butyl	7"	9"	.015"
800-GB/9/10	One-Piece Glove	Butyl	10"	9"	.015"
800-GH	One-Piece Glove	Hypalon	8"	9"	.015"
800-GH/7	Two-Piece Glove	Hypalon	8"	7"	.015"
800-GH/8	One-Piece Glove	Hypalon	8"	8"	.015"
800-GH/10	One-Piece Glove	Hypalon	8"	10"	.015"
800-GH/8/10	One-Piece Glove	Hypalon	10"	8"	.015"
800-GH/9/10	One-Piece Glove	Hypalon	10"	9"	.015"
800-GH/9/6	One-Piece Glove	Hypalon	7"	9"	.015"
<b>815-PGB, "La Petite" Glove Box Gloves</b>					
800-GLLP	One-Piece Glove	Latex	7"	9"	.015"
800-GLLP/7	One-Piece Glove	Latex	7"	7"	.015"
800-GLLP/10	One-Piece Glove	Latex	7"	10"	.015"
800-GN	One-Piece Glove	Neoprene	8"	9"	.015"
800-GN/8	One-Piece Glove	Neoprene	8"	8"	.015"
800-GN/10	One-Piece Glove	Neoprene	8"	10"	.015"
800-GN/11	One-Piece Glove	Neoprene	8"	11"	.015"
800-GN/9/6	One-Piece Glove	Neoprene	7"	9"	.015"