



ID 160

Injector nozzle
with plastic supports, Ø 4, length 160 mm.

- Suitable for standard mobile injector units and injector baskets
- Suitable for measuring flasks and test tubes, for example



Product details - ID 160

EAN: 4002511116993 / Article no. 69716001D / Mat. no. 03810350

Product category		^
Washer-disinfectors, laboratory		•
Large-capacity washer-disinfectors, laboratory		•
Machine/accessory categorisation		^
G 7825		•
PG 8527		•
PG 8535		•
PG 8536		•
PG 8583		•
PG 8583 CD		•
PG 8593		•
Product category		^
Injector nozzle		•
Product properties		^
Material	Stainless steel Plastic	
Colour	Stainless steel Black	
Area of application		^
Reprocessing laboratory glassware		•
Dimensions and weight		^
External dimensions, depth in mm	160	
Diameter in mm	4	

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:



Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432
Fax: (01) 4523967
E-mail: info@labunlimited.com
Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30
Fax: 08452 30 50 30
E-mail: info@labunlimited.co.uk
Web: www.labunlimited.co.uk