

LABORHITS

LaboPlast® Laboratory stands, PP



Easy and flexible height adjustment

Comfortable and smooth operation

With innovative magnet concept

- White, shatter-proof, resistant to chemicals
- For 2 or 4 funnels. Sturdy and stable funnel stand
- Adjustable height up to 450 mm
- Thanks to an innovative magnet concept no clamping screws are necessary. The height of the funnel holder can be adjusted easily and flexibly with one hand
- With 2 or 4 reducing inserts for small funnels
- For funnels von Ø 40 - 180 mm

	Funnel stands	width mm	length mm	height mm	PK	Cat. No.
1	2 funnels	138	340	500	1	9.069 922
2	4 funnels	251	251	500	1	9.069 924

LABORHITS

Stand for Imhoff-sedimentation funnel, PP



The stand for sedimentation funnels is used to keep, hold or store Imhoff funnels safely and reliably during sedimentation. Two funnel holders ensure that Imhoff funnels will be held safely and securely and can be placed in an absolutely vertical position. Stand with two holders for sedimentation funnels, suitable for funnels 90 mm or more in diameter. Thanks to an innovative magnet concept no clamping screws are necessary. The height of the funnel holder can be adjusted easily and flexibly with one hand. With the corresponding accessory the holders for Imhoff funnels can also be replaced by holders for all-purpose funnels.

- PP, chrome-plated steel
- Adjustable height up to 450 mm
- For funnels with Ø 90 mm or more



Description	PK	Cat. No.
Stand for 2 sedimentation funnels	1	6.201 027

Accessories for funnel stand LaboPlast®, PP

Description	For	PK	Cat. No.
1 Funnel holder with reducing insert	Stand with 2 or 4 funnels/ Imhoff funnels	1	6.284 549
2 Sedimentation funnel holder	Stand with 2 or 4 funnels/ Imhoff funnels	1	6.284 550
3 Replacement set for 40 pipettes	for stand with 4 funnels	1	6.284 551