
Data sheet

Neomycin sulfate

Cat. No: AB018 Cat. No: AB019 Cat. No: AB021 Cat. No: AB022

Introduction

Neomycin sulfate (NS) is an aminoglycoside antibiotic that is active against Gram-negative and many Gram-positive bacteria as well as some protozoa. NS inhibits translation by binding to the small subunit of prokaryotic ribosomes and is a calcium channel protein inhibitor.



Product Size

	AB018	AB019
Neomycin sulfate(powder)	1g	5g
	AB021	AB022

Storage conditions

Store at -20° C

Application

- ✓ Cell culture
- Plasmid selection
- ✓ Gene regulation

Technical specification

- Formula: C23H46N6O13.xH2SO4
- CAS: 1405-10-3
- Molecular weight: 614.7 (Base)
- Appearance: White to off-white solid.
- Potency: min. 680
- pH (1% in H₂O): 5.0 7.5
- Sulfate (%): 27.0 31.0
- Sulfated ash: <1.0
- Solubility: Soluble in water

PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively for research purposes and in vitro use only. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals. Please refer to www.canvaxbiotech.com for Material Safety Data Sheet of the product.



Lab Un

Tallaght Business Park Whitestown, Dublin 24, Ireland D24 RFK3

Tel: (01) 4523432 Fax: (01) 4523967 Web: www.labunlimited.com

Quatro House, Frimley Road, Camberley, United Kingdom

E-mail: info@labunlimited.com E-mail: info@labunlimited.co.uk Web: www.labunlimited.co.uk

GU16 7ER Tel: 08452 30 40 30 Fax: 08452 30 50 30





Notes

Distributed by:



Tallaght Business Park Whitestown, Dublin 24, Ireland D24 RFK3

Tel: (01) 4523432 Fax: (01) 4523967 Web: www.labunlimited.com

Quatro House, Frimley Road, Camberley, United Kingdom GU16 7ER

Tel: 08452 30 40 30 Fax: 08452 30 50 30 E-mail: info@labunlimited.com E-mail: info@labunlimited.co.uk Web: www.labunlimited.co.uk

