

**O 184**

Upper basket
for optimum loading of centrifuge tubes, phials
or test tubes, for example.

- Equipped with 96 injector nozzles
- Equipped with connection for hot air drying
- Vertical clearance: 170 mm
- Powder dispensing not possible



Product details - O 184

EAN: 4002511114937 / Article no. 69118401D / Mat. no. 03807480

Product category ^	
Washer-disinfectors, laboratory	•
Machine/accessory categorisation ^	
PG 8535	•
PG 8536	•
Product category ^	
Basket for vials	•
Basket for test tubes	•
Upper basket	•
Product properties ^	
Material	Stainless steel Plastic
Colour	Stainless steel Black
Maximum vertical clearance in mm	170
Area of application ^	
Reprocessing pipettes/vials	•
Capacity ^	
Test tubes [number]	96
Dimensions and weight ^	
External dimensions, height in mm	468
External dimensions, width in mm	475
External dimensions, depth in mm	531

Net weight in kg	4,500
------------------	-------

Standard accessories	^
----------------------	---

ID 90, injector nozzles Ø 2.5 mm, length 90 mm [number]	96
---	----

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:

Lab Unlimited
CARL STUART GROUP

Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432
Fax: (01) 4523967

E-mail: info@labunlimited.com
Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30
Fax: 08452 30 50 30

E-mail: info@labunlimited.co.uk
Web: www.labunlimited.co.uk