



U 175/1

Lower basket
for the optimum loading of narrow necked
glassware.

- Equipped with 33 injector nozzles
- Vertical clearance: 170 mm¹
- Only suitable for use with upper baskets O 175 / O 187 / O 184



Product details - U 175/1

EAN: 4002513805543 / Article no. 69117503D / Mat. no. 06439170

Product category		^
Washer-disinfectors, laboratory		•
Machine/accessory categorisation		^
PG 8536		•
Product category		^
Lower basket		•
Product properties		^
Material	Stainless steel	
Colour	Stainless steel	
Maximum vertical clearance in mm	170	
Area of application		^
Reprocessing laboratory glassware		•
Capacity		^
Narrow-necked glassware [number]	33	
Dimensions and weight		^
External dimensions, height in mm	250	
External dimensions, width in mm	475	
External dimensions, depth in mm	531	
Net weight in kg	6,400	
Standard accessories		^
F 351 injector nozzles Ø 4 mm, length 160 mm [number]	33	

E 353 clips for injector nozzles E 351 [number]

33

Rinse nozzle for powder dispenser [number]

1

Subject to technical changes; no liability accepted for the accuracy of the information given!

¹ Depending on the laboratory glassware the max. vertical clearance can be lower.

Distributed by:

Lab Unlimited
CARL STUART GROUP

Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432
Fax: (01) 4523967
E-mail: info@labunlimited.com
Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30
Fax: 08452 30 50 30
E-mail: info@labunlimited.co.uk
Web: www.labunlimited.co.uk