



Y S I Environmental

Distributed by:

Lab Unlimited
CARL STUART GROUP

Tallaght Business Park
Whitestown
Dublin 24
Ireland

Tel: (01) 4523432
Fax: (01) 4523967
E-mail: info@labunlimited.com
Web: www.labunlimited.com

Unit 59, Frimley House
35 High Street
Frimley, Surrey
GU16 7JQ

Tel: 08452 30 40 30
Fax: 08452 30 50 30
E-mail: info@labunlimited.co.uk
Web: www.labunlimited.co.uk

YSI BOD Analyst Pro™

Desktop Software for the calculation of BODs - Easy, Efficient, Accurate

- Automatically calculates BODs/CBODs conveniently and easily
- Improve data integrity by eliminating human error in calculations
- Save time by quickly creating daily bench sheets from established samples and groups
- Use with Model 52, 5000, 5100, Professional Plus or ProODO® instruments for direct connectivity; manual entry also possible
- All batches are saved for simple retrieval at any time
- User-configurable set up regardless of seeding method, number of dilutions, or other test criteria
- Bar code compatibility ensures efficiency and improves accuracy

+1 937 767 7241
800 897 4151 (US)
www.yssi.com
environmental@yssi.com

YSI Gulf Coast
+1 225 753 2650
environmental@yssi.com

YSI Hydrodata (UK)
+44 1462 673 581
europe@yssi.com

YSI Middle East (Bahrain)
+973 3977 1055
halsalem@yssi.com

YSI (Hong Kong) Limited
+852 2891 8154
hongkong@yssi.com

YSI (China) Limited
+86 10 5203 9675
beijing@yssi-china.com

YSI Nanotech (Japan)
+81 44 222 0009
nanotech@yssi.com

YSI South-Asia
+91 989 122 0639
sham@yssi.com

YSI Australia
+61 7 390 17223
acorbett@yssi.com

YSI Brazil
+55 48 9942 7840
psterling@yssi.com

Who's Minding the Planet?, We Know D.O., ProODO and BOD Analyst Pro are registered trademarks of YSI Incorporated.

©2010 YSI Incorporated
Printed in USA W13 0910

Specifications are subject to change. Please visit www.YSI.com to verify all specs.

ISO 9001
ISO 14001



YSI incorporated
Who's Minding the Planet?®



The main screen is very simple to navigate. Quickly set up samples or groups and instantly build batches for the days BODs.

YSI BOD Analyst Pro Specifications

Instrument	Model 52, 5000, 5100, Professional Plus or ProODO direct connectivity;
Compatibility	Initial and final values can also be manually entered from any DO instrument; can be used without direct connectivity
Operating System	Microsoft Windows® XP, Vista or 7

YSI BOD Analyst Pro Ordering Information

Item #625120	Includes CD with User Manual; compatible with YSI instruments 52, 5000, 5100, Professional Plus and ProODO optical dissolved oxygen for direct instrument to PC communication
--------------	---

YSI.com/BOD