



Specifications	Lab Unlimited GF 1000 Kit
Order info:	Z-CSL-100-597
Mode of Measurement	Displays mass flow, volumetric flow, temperature, and pressure readings simultaneously
Split Ration Mode	N/A
Gases	Various gas options, 30+
Flow Range	0-1000 SCCM and mL/min (other ranges available upon request)
Resolution	0.1 mL/min and 0.1 SCCM
Accuracy	±(0.8% of reading +0.2% of full scale)
Display	TFT-LCD
Calibration	NIST traceability
Operating Temperature Range	-10°C to 60°C
Size	60.33 x 170.43 x 23.5 mm
Weight	453.59g (1lb)
Power	Rechargeable battery with universal USB micro charger
CE mark	Yes
Data output	RS232 Output
Auto Power off	N/A
Inlet ID/OD	1/8" Female NPT thread
Inlet Pressure	Max 145psig recommended 50psig
Additional Features	Turndown Ratio of 1000:1 Measures pressure (inlet and Vent of GC)