

PT31 | Peltier Thermostat

This electronic water-bath thermostat with Peltier element is a versatile, high-performance instrument.

In one application, for example, it can be used to set the correct refractometer temperature. It is extremely robust, compact and easy to operate. Because it is so small it does not take up valuable space in the laboratory.

Resolution	0.1 °C
Heating power	30 W
Cooling power	15 W
Power supply	115–230 V
Pump pressure	2000 Pa
Pump performance	20 l/h
Temperature	8–40 °C (continuously adjustable)
Temperature accuracy	±0.2 °C
Bath volume	ca. 100 ml
Dimensions in cm	8.0 x 21.0 x 14.0
Weight	1.5 kg



Further Products available from A.KRÜSS Optronic GmbH

- Refractometers
- Polarimeters
- Density Meters
- Microscopes
- Melting Point Meters
- Cold Light Sources
- Spectroscopes
- Spectral Measuring Instruments
- UV Lamps
- Software
- Gemmological Instruments

Distributed by:

Lab Unlimited
CARL STUART GROUP

Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432
Fax: (01) 4523967
E-mail: info@labunlimited.com
Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30
Fax: 08452 30 50 30
E-mail: info@labunlimited.co.uk
Web: www.labunlimited.co.uk