

# Multichannel pipette ME and MT



## Multichannel pipette ME/MT

The fully autoclavable PIPETTES ME/MT are specially designed for working with gloves under sterile conditions.

The volume range from 1µl to 300 µl is covered by 4 models each in the 8- and 12-channel versions.

The patented manifold system ensures individual tight sealing of the pipette tips on each channel.

- perfect measuring results
- fully autoclavable and UV-light resistant
- volume settings by pushbutton avoids snatching the gloves
- ergonomic design, handle can be turned to an individual convenient position
- patented system makes manual corrections of the tip sealing redundant
- maintenance free: all components are manufactured either from stainless steel or chemical resistant plastics
- each pipette comes with quality certificate, holder and calibration key

### 8-channel-pipette ME

Volume range µl	Model	Accuracy max. vol %	Precision < %	Quantity per Pack	Order No.
1-10	ME10	2,0	1,2	1	47 70 064
5-50	ME50	1,6	0,6	1	47 70 065
20-200	ME200	1,0	0,6	1	47 70 066
50-300	ME300	1,0	0,6	1	47 70 067

### 12-channel-pipette MT

Volume range µl	Model	Accuracy max. vol %	Precision < %	Quantity per Pack	Order No.
1-10	MT10	2,0	1,2	1	47 70 068
5-50	MT50	1,6	0,6	1	47 70 069
20-200	MT200	1,0	0,6	1	47 70 070
50-300	MT300	1,0	0,6	1	47 70 071

Distributed by:

**Lab Unlimited**  
CARL STUART GROUP

Tallaght Business Park  
Whitestown, Dublin 24,  
Ireland  
D24 RFK3

Tel: (01) 4523432  
Fax: (01) 4523967  
E-mail: [info@labunlimited.com](mailto:info@labunlimited.com)  
Web: [www.labunlimited.com](http://www.labunlimited.com)

Quatro House, Frimley Road,  
Camberley,  
United Kingdom  
GU16 7ER

Tel: 08452 30 40 30  
Fax: 08452 30 50 30  
E-mail: [info@labunlimited.co.uk](mailto:info@labunlimited.co.uk)  
Web: [www.labunlimited.co.uk](http://www.labunlimited.co.uk)