

661 Route 23 South
Wayne, New Jersey 07470
Tel: 973-694-0500 \* Fax: 973-694-7199
Internet Address: http://www.bel-art.com

Date: 2/22/2017

Re: Expiration Dates / Sterilization

**Sterileware® Product Line** 

Bel-Art products do not carry expiration dates. The products remain sterile as long as the barrier remains unbroken preventing bacteria and or other airborne particulates from entering into the sterile item.

It is recommended that the items be properly stored to maintain the packages integrity. Items should be stored in a clean and dry (non-humid) environment while protected from direct sunlight until ready for use. It is also recommended that the principal processing technique of first in – first out be used while maintaining inventory. Any definitive shelf life should be governed by the end-users protocol.

Devices are generally regarded as sterile as long as package integrity has not been compromised.

I certify that I am authorized to declare that the statement made with this Certificate of Compliance is true and accurate.

Signature:

Michael Grzelak

SP Scienceware, a Divisions of SP

661 Route 23 South, Wayne, NJ 07470

O: +1.973.694.0500 x 4217 mgrzelak@belart.com